

The specialists' first choice



SG3 blade & No 3 Stainless handle

A comprehensive specialist range of Skin Graft blades and handles, NOW including blades for Dermatomes

The Swann-Morton range now offers the knives required by a surgeon to harvest both "split-thickness" and "full-thickness" skin grafts used in all areas of reconstructive surgery.

From a standard scalpel to Skin Graft Knives all manufactured to the highest quality standards and with the precision and consistency of performance long associated with the Swann-Morton name.

The SG3 "Minor" Sterile Skin Graft blades can be used conveniently with a standard No.3 Surgical handle at the site of a major incident where time may be critical and debridement of the wound could result in improved healing times and

reduction in infection of the surrounding healthy tissues.

This range provides the surgeon with a selection of three Stainless Steel Skin Graft Handles, namely the Braithwaite, Cobbett and Watson that can be used with our Sterile 158mm Stainless Steel Skin Graft blades on larger donor sites.

Plus NEW to the range are our choice of 3 Stainless Steel Sterile Dermatome blades which are used primarily within Plastic, Reconstructive and Oral surgery procedures making Swann-Morton the "Specialists" first choice.

Swann-Morton®
Made in Sheffield - Precisely



CE

All medical devices are CE marked in accordance with the Medical Device Directive (93/42/EEC)



FM 73368



Made in Sheffield - *Precisely*

THE WORLD'S LEADING SURGICAL BLADE MANUFACTURER

The 'NEW' range of Dermatome blades

Dermatome blades - Sterile Stainless Steel			Product	Product No	Unit
Product	Product No	Unit	2 Segmental	9944	Box of 10
1 Segment	9943	Box of 10	3 Segmental	9945	Box of 10

These blades are individually wrapped in foil packets and boxed in units of 10.



Skin Graft handles, blades and sets



Skin graft blades - Sterile Stainless Steel.			Braithwaite skin graft knife set		
Product	Product No	Unit	For use with Blade 9940		
SG3*	2201	20 blades	Product	Product No	Unit
158mm	9940	10 blades	Braithwaite Skin Graft Knife Set	9904	1 set
Braithwaite skin graft knife - Non-sterile Stainless Steel, Non-Sterile, individually boxed.			Cobbett skin graft knife set		
For use with Blade 9940			For use with Blade 9940		
Product	Product No	Unit	Product	Product No	Unit
Braithwaite Skin Graft Knife	9901	1 handle	Cobbett Skin Graft Knife Set	9905	1 set
Cobbett skin graft knife - Non-sterile Stainless Steel, Non-Sterile, individually boxed.			Watson skin graft knife set		
For use with Blade 9940			For use with Blade 9940		
Product	Product No	Unit	Product	Product No	Unit
Cobbett Skin Graft Knife	9902	1 handle	Watson Skin Graft Knife Set	9906	1 set
Watson skin graft knife - Non-sterile Stainless Steel, Non-Sterile, individually boxed.			All sets contain Stainless Steel Handle with 2 boxes of 10 Sterile 158mm Stainless Steel skin graft blades.		
For use with Blade 9940			*Fits standard No. 3 handle		
Product	Product No	Unit	For complete product specifications and comprehensive cleaning instructions, go to www.swann-morton.com		
Watson Skin Graft Knife	9903	1 handle			

Swann-Morton and the 'ring pattern' logo are registered trademarks of Swann-Morton Limited and related companies. The products shown are manufactured for use by trained health care professionals and can cause injury if used or handled incorrectly. The information and illustrations used in connection with the products shown on this leaflet are for general guidance only and users must seek proper training and instruction in their use before using or handling such products. Swann-Morton Limited exclude all liability for injury to the extent permitted by law.



Swann-Morton Ltd. Penn Works, Owlerton Green, Sheffield S6 2BJ
 Telephone: 0114 2344231, Sales 0114 2344223
uksales@swann-morton.com, exportsales@swann-morton.com
www.swann-morton.com